



01 WATER PEELING

02 ION TREATMENT

03 LIGHT THERAPY

04 GOLD THERAPY

05 ION ACTIVATION

REJUVENA™
PREMIUM

5 SECRET TIPS TO MAINTAIN YOUR BEAUTY

REJUVENA™ PREMIUM



RICHWAY & FUJI BIO INC. 1314 SOUTH KING STREET SUITE 520, HONOLULU, HI 96814 USA
TEL: (808) 589-2800 TOLL-FREE: (855) 338-6410 FAX: (808) 597-1651
WWW.RICHWAYANDFUJIBIO.COM

SHAOXING RICHWAY TRADING CO., LTD. ROOM 2703, ZHONGYANG BUILDING, YUMIN ROAD,
KEQIAO DISTRICT, SHAOXING, ZHEJIANG, CHINA 312000
TEL: +86-575-8110-1350 FAX: +86-575-8110-1370

RICHWAY EUROPE B.V. BEECHAVENUE 54-80, 1119 PW SCHIPHOL-RIJK, NETHERLANDS
TEL: +31-(0)20-658-6231/+31-(0)20-303-3411 (INTERNATIONAL USE) FAX: +31-(0)20-658-6142

RICHWAY & LIFE JAPAN CORP. 2-6-41 NAKAHATA, ONOJO-SHI, FUKUOKA,
JAPAN 816-0921 TEL: +81-92-575-1004 FAX: +81-92-575-1006

R&L CORPORATION 11F, ACE GWANGMYEONG TOWER B, 108 HAAN-RO,
GWANGMYEONG-SI, GYEONGGI-DO, KOREA 14319
TEL: +82-2-6112-7700 FAX: +82-2-6112-7749 WWW.RICHWAYKOREA.COM

PUBLISHED 09/2016

REJUVENA™ PREMIUM

PEELING
ION TREATMENT
LIGHT THERAPY
3-IN-1



THE
RICHWAY
GROUP

RICHWAY & FUJI BIO INC.
1314 South King Street Suite 520
Honolulu, HI 96814 USA
Tel: 808-589-2800 Fax: 808-597-1651
www.richwayandfujibio.com

The Secret Behind Beautiful Skin

REJUVENA™ PREMIUM

Use with your daily basic skincare products for 15 minutes a day to remove oil, dirt, and debris from pores and to deeply nourish your skin with a gold massage and ion treatment. Richway & Fuji Bio has developed the Rejuvena™ Premium to respond to the real needs of aging skin. You will feel as if you were just treated by a professional aesthetician.



REFRESH AND REJUVENATE
A SPA-LIKE TREATMENT
AT HOME!



AESTHETIC TREATMENT
IN THE PALM OF YOUR HAND

3-STEP SYSTEM

PEELING : MOISTURIZING + PEELING

G-BLUE : NOURISHING + BLUE LIGHT THERAPY

G-RED : NOURISHING + RED LIGHT THERAPY

REJUVENA™ PREMIUM

5 Main Points



Water Peeling

Gentle peeling and oil and dirt removal

Gentle vibrations will cleanse your pores of dirt and debris with water. This step helps remove dead skin cells to make your skin both clear and clean.



Ion Treatment

Nutrient absorbency and sonic massage

Ions penetrate nutrients deep into the skin and microcurrent treatment on the surface of the skin gives you a luminous glow.



Light Therapy

Skin-soothing effects and improved elasticity

Blue (470 nm) & Red (624 nm) Light treatment refreshes and revitalizes exhausted skin.

BLUE: Reduce visible signs of aging and sagging skin.

RED: Skin rejuvenation treatment.



Gold Therapy

Skin beautifying and youthful glow

Gold massage maximizes light-therapy results to achieve glossy skin.



Ion Activation

Positive and negative ions

The ions will be activated by placing your hand along the gold contact strips.

SIMPLE AND EASY

One-Button System

When you press and hold the power button while the device is in OFF mode, the device will turn on with a sound to alert you that the device is now ON. When you press and hold the power button while the device is in use or in standby mode, the device turns off with a sound.

*** The device switches to standby mode automatically after 5 minutes of use. The device turns off automatically after 5 minutes of inactivity in standby mode.**



STEP 1. GENTLE PEELING

After washing, lightly touch wet skin with the blade to remove impurities and dead skin cells.



STEP 2. G-BLUE (470 nm)

Release negative ions into the skin along with a sonic massage after applying a skin ampoule or an essence to deeply infuse nutrients from the skincare products. This alleviates skin troubles.



STEP 3. G-RED (624 nm)

Release positive ions into the skin along with a sonic massage after applying nutrients or moisturizing cream to maintain skin elasticity.



Low Battery and Charging

If the product runs out of battery power while in use, 3 LED lights will blink 5 times and then the product will shut off automatically. Please connect the Micro USB Charger into the charging socket located on the bottom of the device. Once the device's battery begins to charge, the red LED light will turn on. The blue LED light will turn on once the product has finished charging. After the product has been fully charged, please remove the charging cable before continuing to use the product.